



ConfigFile explanation (1 dataset = 1 Config file)

Overview : 1 dataset = 1 Config file

This page will help you to understand and set all the values of the config file

`DeepCellMap/code/config/dataset_new_fluorescence.py` or

`DeepCellMap/code/config/dataset_new_ihc.py` depending on the data on hand.

Pour moi :

Je vais voir les attributs les uns apres les autres

Je les map dans le following table et je copie cole le "how to choose" en commentaire avant l'attribut du config file

☐ To add to the other config files

☐ mapping_cells_colors (deja dans ihc)

☐ limit_row_col_for_slide (deja dans ihc)

Common attributes

Attribute name	How to choose	IHC	Fluorescence
dataset_name	Use a descriptive name for your dataset.	✓	✓
mapping_img_name	Specify the names of the images in the dataset.	✓	✓
mapping_img_gender	Specify the genders of the individuals in the dataset.	✓	✓
mapping_img_disease	Specify the disease of the individuals in the dataset.	✓	✓
mapping_img_age	Specify the ages of the individuals in the dataset.	✓	✓
cell_class_names	List of cell/object names	✓	✓
mapping_cells_colors	Color associated to each cell	✓	✓
debug_mode	Set to <code>True</code> if you want to enable debugging mode.	✓	✓
data_type	Specify the type of data, such as whole slide images (WSI) or fluorescence	✓	✓
consider_image_with_channels	A fluorescent image can be considered sometimes as RGB (eg Cancer data) then,	✓	✓
data_format	Specify the file format of the data, such as TIFF.	✓	✓
conversion_px_micro_meter	Specify the conversion factor from pixels to micrometers. Ex : 1 px = 0.45 micro meter	✓	✓
channel_names	Specify the names of the channels in the data.	✗	✓
dim_position	In tensor of the image, position of C, X, Y and Z	✗	✓
has_Z	If there is a Z canal	✗	✓

tile_width	Set the width of each image tile.	✓	✓
tile_height	Set the height of each image tile.	✓	✓
roi_border_size	Set the size of the border around the region of interest (ROI).	✓	✓
border_size_during_segmentation	Set the size of the border during segmentation.	✓	✓
crop_size	Set the size of the crop for each image tile.	✓	✓
scale_factor	Downscaling factor	✓	✓
threshold_tissue	Set the threshold for tissue segmentation. If a tile contains at least (1-threshold_tissue)% of tissue the tile is considered to have tissue inside. (Save computation time)	✓	✓
preprocessing_config	Configure the preprocessing parameters for the data.	✓	✓
tissue_segmentation_param	Set the parameters for tissue segmentation.	✓	✓
tissue_extraction_accept_holes	If the segmented tissue can contains holes	✓	✓
save_tissue_segmentation_steps	Specify whether to save the intermediate steps of tissue segmentation.	✓	✓
channel_used_to_segment_tissue	Which channel is used to segment the tissue	✗	✓
use_imgs_as_channels	Some images can be split in different images (channel 1 is 1 image, channel 2 is one image and so on). Specify whether to use images as channels. So that image recombination is done. Use for Cancer dataset.	✓	✓
channels_cells_to_segment	List of the channels to segment during the process of cells detection (ex [1,2])	✗	✓
channels_of_interest	List channels of interest during the application of DeepCellMap	✗	✓
cells_from_multiple_channels	Only fluorecence : usually 1 channel = 1 celltype but cells can also result from the combination of several channels. Ex : dict("iba1_cd68" : [1,2]) → new celltype iba1_cd68 emerges from the combination of channel 1 and 2. Cells that are double positive are considered iba1_cd68	✗	✓
association_cell_name_channel_number	dict specifying the association between cell names and channel number. Cells created from several channels are associated the channel number max_channel + key_position in cells_from_multiple_channels	✗	✓
param_best_cellpose	Set the parameters for the best Cellpose model (to detect ovoid nuclei)	✓	✓
cellpose_parameters	Dict that contains main Cellpose parameter like <code>downscale_factor_cellpose_tissue_segmentation</code> , <code>channel_nuclei</code> and <code>tile_subdivision_factor</code> used during the creation of the density heatmap	✓	✓
model_segmentation_name	Set the name of the segmentation model.	✓	✓
cell_segmentation_param	Set the parameters for cell segmentation.	✓	
cell_segmentation_param_by_channel	Set the parameters for cell segmentation but by channels (and can be image-specific) first key is		✓

	default parameters and other keys are the parameters for the other images		
classification_param	Set the parameters for cell classification.	✓	✗
cell_class_names_for_classification	Name of all the categories used during classification (include Background, and Detected)	✓	✗
tile_test_segmentation	Set the parameters for tile-level segmentation.	✓	✓
roi_test_tissue_border	Set the parameters for testing tissue border within ROIs.	✓	✓
roi_cool	Set the parameters for cool ROIs.	✓	✓
statistics_with_proba	Specify whether to compute statistics with probabilities. Used when DL did the classification and gave each cells a probability to belong to each celltype	✓	✗
colnames_table_cells_base	Specify the column names for the table of cells (created during detection segmentation and classification)	✓	✓
colnames_df_image	Specify the column names for the image dataframe.	✓	✓
colnames_df_roi	Specify the column names for the ROI dataframe.	✓	✓
path_cells	Set the path to the directory where cell images are stored.	✓	✓
physiological_regions_max_square_size	Set the maximum square size for each physiological region.	✓	✗
physiological_regions_group_for_comparaison	Set the group for comparison for each physiological region.	✓	✗
cell_cell_colocalisation_config	Set the parameters for cell-cell colocalisation analysis.	✓	✓
dbscan_based_analysis_config	Set the parameters for DBSCAN-based analysis.	✓	✓
neighbors_analysis_config	Set the parameters for neighbors-based analysis.	✓	✓
limit_row_col_for_slide	Limite row and col (on some slides the tissue is repeated)	✓	✓